




3M Food Safety

Testing Supply Technologies for Environmental Monitoring

Integrated environmental monitoring

Today you'll learn about:

- Importance of environmental monitoring
- Objectives of environmental sample collection and a risk-based approach to site selection
- Guidance for building components of a holistic program, including testing for ATP, indicator and spoilage organisms, pathogens, and allergens
- Best practices for data trending



3M.com/EnvironmentalMonitoring




How to focus on preventing food safety issues

Reactive approach, waiting for final product status to trigger action.



Proactive approach, establishing interventions to minimize risk of contamination.



Shifting testing priorities from finished product testing to a more strategic approach that identifies contamination sources.

Controlling distribution systems to maintain the food safety and quality.

What is environmental monitoring?

Sampling and testing the environment and equipment within a food manufacturing facility to prevent cross contamination of the finished product from the environment.

Environmental monitoring should be approached holistically, encompassing a range of tests to ensure both food safety and quality.



ATP Monitoring



Indicator
Organisms



Pathogens



Spoilage
Organisms



Allergens

Why focus on environmental monitoring?

Food safety events

Food safety events are often attributed to failures of Prerequisite Programs (PRPs) and not Critical Control Points (CCPs).

- Sanitation
- Sanitary design
- Personnel Hygiene
- Raw material storage

Regulations and standards


Regulations and standards are including validation and verification of activities previously considered Prerequisite Programs.

- FDA FSMA requires environmental monitoring as a sanitation verification activity, if post-lethality contamination is a risk

Financial return on investment

- Reduce risk of product recalls
- Limit the size and scope of recalls
- Mitigate brand damage
- Possibly validate extended production run times
- Increase overall productivity in the plant

Guidance on building an environmental monitoring program



Identify sampling sites

Determine sampling frequency

Analyze the data


Perform root cause analysis

Implement Corrective Actions

Selecting sampling sites using zones

	Zone 1 Food Contact Surfaces Slicers, peelers, fillers, hoppers, screens, conveyor belts, air blowers, employee hands, knives, racks, work tables
	Zone 2 Non-Food Contact Surfaces in Close Proximity to Food and Food Contact Surfaces Processing equipment exterior and framework, refrigeration units, equipment control panels, switches
	Zone 3 More Remote Non-Food Contact Surfaces Located In or Near the Processing Areas Forklifts, hand trucks, carts, wheels, air return covers, hoses, walls, floors, drains
	Zone 4 Non-Food Contact Surfaces Outside of the Processing Areas Locker rooms, cafeterias, entry/access ways, loading bays, finished product storage areas, maintenance areas

———— Higher focus/frequency
- - - - Lower focus/frequency




Test point selection process:

How do you determine which sites to swab?

Risk-Based Approach

How hard is it to clean?

How close to the food?



		Probability		
		Low	Medium	High
Hazard	High	Once per week/month	High frequency testing	High frequency testing
	Medium	Low frequency testing	Once per week/month	High frequency testing
	Low	Low frequency testing	Low frequency testing	Once per week/month

Integrated environmental monitoring.


Not just ATP...

Not just Pathogens...

Not just Allergens...

Not just Indicator Organisms...

A holistic EMP should address
relevant hazards in a plant!



Sample to find

Growth Niches

- Location that supports microbiological growth and is protected from the sanitation process
- Characterized by high microbial counts after cleaning and sanitation
- May trap nutrients or water which can promote the growth of bacteria
- Can be a potential place for the formation of biofilms



Selection of the sample collection device

Select a sample collection device that:

Reaches the sample target area

Aseptically collects the sample

Dislodges microorganisms from the surface

Neutralizes residual sanitizer

Selection of the sample collection device



Sponges

Larger sampling devices and are available in a variety of formats.

Preferred choice when the area being sampled is large and readily accessible (greater than 100 cm²).

Preferred if qualitative pathogen testing is to be conducted.



Swabs

Smaller sampling devices consisting of a tip or bud for collecting the sample attached to a long flexible stem.

Due to their smaller size and ease-of-use for sampling a defined area, they are preferred for small crevices and penetrations (areas of 100 cm² or less).

Can be particularly useful for quantitative environmental testing (e.g. for indicator organisms).

Integrated environmental monitoring


Not just ATP...

Not just Pathogens...

Not just Allergens...

Not just Indicator Organisms...

A holistic EMP addresses relevant hazards in a plant!



The importance of hygiene monitoring



Cleaning and sanitation to remove or prevent pathogens, spoilage microorganisms and allergens.

Hygiene monitoring to assess the effectiveness of your cleaning and sanitation programs.

Just because it looks clean,
doesn't mean it is clean

Microbiological
testing
(swab & plate)

Visual Inspection


ATP
bioluminescence

All have their place in a robust hygiene monitoring program

Hygiene monitoring generates data to determine if a surface has been cleaned sufficiently
and enables the high-risk decision to start food production

ATP as a tool for hygiene monitoring

Principle behind ATP bioluminescence



ATP = Adenosine Triphosphate
The “energy currency” molecule
of all living organisms




Luciferin/luciferase



Increase
In organisms or
organic residues

Increase
In ATP levels

Increase
In light (RLU)



Sample site selection example: meat slicer

Guard			
	Low	Medium	High
High	★		
Medium			
Low			

Blade			
	Low	Medium	High
High			★
Medium			
Low			

Back Plate			
	Low	Medium	High
High		★	
Medium			
Low			

Slicer Handle			
	Low	Medium	High
High			
Medium	★		
Low			

Collection Area			
	Low	Medium	High
High	★		
Medium			
Low			

Slice Thickness Knob			
	Low	Medium	High
High			
Medium		★	
Low			



Integrated environmental monitoring


Not just ATP...

Not just Pathogens...

Not just Allergens...

Not just Indicator Organisms...

A holistic EMP addresses relevant hazards in a plant!



“I use ATP, why would I need to do indicator testing?”

ATP verifies
Cleaning

Indicators validate and verify
Sanitization

Organisms can exist in sites where sanitizer is not able to reach.

If the processing environment is under control, indicator organisms will be under control.



Spoilage and indicator organisms

Indicators

Help us to determine hygienic status of the processing equipment and environment.

Validate and verify sanitation and process control steps.

- Total Plate Count
 - Coliforms
 - *Enterobacteriaceae*
-

Spoilage

EMP allow companies to take a proactive approach to microbial spoilage, rather than retrospectively addressing failures as they arise.

- Total Plate Count
 - Yeast and Molds
 - Lactic Acid Bacteria
-


Identify sources of contamination and sampling sites

Sources

- Cross-contamination of sources outside the production area or in process (water)
- Ingredients and raw materials
- Atypical activities (change in cleaning chemicals, new personnel)
- Equipment design/repairs
- Part wear

Frequency and application

- After each sanitization cycle and before starting production (useful strategy to perform data trending)
- Pre-production: run the equipment (conveyor belts) before sampling to have a better chance of collecting residual microorganisms that remained after the sanitization process
- During production: ideally 3-4 hours into run or at multiple time points
- After non-routine activities such as construction or maintenance




Continuous improvement

Monitoring for indicator and spoilage organisms is beneficial in that the quantitative results allow for baseline determination which defines the acceptable limits and analysis of trends.

This allows for:

- Identification of a failure point that triggers an investigation and a root cause analysis
- Action before a failure point is reached
- Understanding of seasonality effects
- Identify opportunities for operational and product improvements

Figure 2. Example of coliform counts and baseline-level post-sanitation



Take the time to analyze results in order to gain the full benefit of implementing environmental monitoring for indicator and spoilage organisms.

Integrated environmental monitoring


Not just ATP...

Not just Pathogens...

Not just Allergens...

Not just Indicator Organisms...

A holistic EMP addresses relevant hazards in a plant!



Key parameters of pathogen environmental monitoring



Risk-Based Approach

What to Test For

What microorganism?
What makes sense based on product/process?
Consider regulatory and customer requirements.

What to Test

Swabs, sponges, dust, air, scrapping material, product residue etc.
Remember sample collection/sample maintenance.

Site Selection

Explore and mapping the entire plant, identify possible sampling sites, potential niches, sites of accumulation, sites with high traffic etc.

Site List

Representation of all possible sampling sites.
It does not imply all have to be sampled every time.

Sample to find

Be aware of the square inch or cm² mentality.

Many training materials and even government guidance documents specify a certain area that should be sampled (typically 12 × 12 in or 30 × 30 cm).

Remember that these are guidelines, as virtually all potential niches that should be sampled as part of an environmental monitoring program are not square or even flat areas.

If sampling sites are not easily accessible, a swab may be more suitable.



If a positive is found

Before start testing for pathogens, you need to get ready for a positive.



Documented

- ✓ Action plan ready for immediate corrective actions by zone.
- ✓ Action plan to verify the absence of a pathogen from an area after a positive finding.
- ✓ Action plan to execute root cause analysis and establish preventive actions.

Recommended to execute
with a presumptive result
(e.g., DNA-based test)

If a positive is found

- Thorough examination of the area
- Visual examination
- Vector swabbing



Use a starburst pattern
Radiate from initial positive
site in all directions (3D).

Collect additional 10–15
samples.

Analyze separately.

Do not composite or pool
the samples.

Integrated environmental monitoring


Not just ATP...

Not just Pathogens...


Not just Allergens...

Not just Indicator Organisms...

A holistic EMP addresses relevant hazards in a plant!



Cross-contact contamination with allergens



Root cause analysis/investigation testing

Zones 1 and 2


How has the cleaning failed?

What needs to be changed to prevent failure again?

Zones 3 and 4

How did the product get there?

People, traffic, construction, spatter, drift, fans, HVAC?



TACT Every day, every time

Time enough time for effective cleaning

Action sufficient mechanical action supplied

Chemical type/concentration

Temperature too hot/too cold

Allergen detection methods

Non-specific test

ATP test swabs



Non-specific protein swabs



- Non-specific
- Indirect measurement
- Application: Cleaning procedure verification of complex (ingredient) production lines

Specific protein test

Quantification and detection of specific protein



Qualitative detection of specific protein



- Highly specific and selective
- Direct measurement of specific proteins
- Application: Validation and verification of cleaning procedures

Integrated environmental monitoring generating knowledge to act



An Environmental Monitoring Program results in an abundant source of data.
What do we do about it?

Data

Facts and figures which rely on something specific,
but which are not organized in any way

120RLU	945RLU	119RLU	120RLU
110RLU	132RLU	987RLU	110RLU
107RLU	156RLU	771RLU	107RLU
605RLU	123RLU	98RLU	605RLU


Integrated environmental monitoring generating knowledge to act




Pattern:
Consistent fails
on weekends

120RLU	945RLU	119RLU	120RLU
110RLU	132RLU	987RLU	110RLU
107RLU	156RLU	771RLU	107RLU
605RLU	123RLU	98RLU	605RLU


Integrated environmental monitoring generating knowledge to act




- Change in staff over the weekend is not following protocol properly
- Additional training or oversight is required



Integrated environmental monitoring generating knowledge to act



- Change in staff over the weekend is not following protocol properly
- Additional training or oversight is required



Trending

Define these to improve clarity.

Who will track and analyze the data?

Who will review the data, and when?

Will upper management support corrective actions?

How will the findings be shared within the company?

Who is accountable, who is a stakeholder?

Consider the big picture of the plant when analyzing data.


Temporary construction project

Changes to activities/sanitation processes/raw material

New equipment/personnel

New QA staff still learning how to sample

Seasonal variation



Supporting a food safety culture

Sharing information enables:

- Staff to see their positive impacts on performance and that goals are being met
- Management to recognize and reward staff
- Management to communicate opportunities for improvement

Supports employee involvement and communication.



Science.
Applied to Life.™




3M Food Safety

Microbiological Solutions and Technology

Agenda


I. Evolution of Quantitative Microbial Methods



II. Evolution of Pathogen Detection Technologies




Microbial enumeration and biochemical differentiation: A very brief timeline.




History of Innovation


3M™ Petrifilm™ Plates




Identify a need: What problem can we solve?




3M Technology Platforms




3M™ Petrifilm™ Aerobic Count Plate construction



A new way to detect gas and confirm results



Coliform is not a taxonomic classification but rather a work definition used to describe a group of Gram-negative, facultative anaerobic rod-shaped bacteria that ferment lactose to produce acid and gas within 48 h at 35°C.



Confirmed Result

Creating the optimal environment

Aerobic


- Products developed for organisms that grow in environments containing oxygen contain films that are gas permeable.

Anaerobic

- Organisms such as lactic acid bacteria grow in the presence of oxygen but thrive in the absence.
- When developing the 3M™ Petrifilm™ Lactic Acid Bacteria Count Plate the inventors sought after new technology to contain a self-contained anaerobic environment.

Anaerobic made easy.

The first self-contained anaerobic plate for lactic acid bacteria testing.



Validations for 3M™ Petrifilm™ Plate Methods



More
than

105

global
validations

International Recognition

AOAC® Performance Tested MethodSM

Environmental *Listeria* Plates
Certificate #030601

Aqua Coliform Count Plates
Certificate #091101

Salmonella Express System
Certificate #061301

Rapid Yeast and Mold Count Plates
Certificate #121301

Rapid Aerobic Count Plates
Certificate #121403

Lactic Acid Bacteria Count Plates
Certificate #041701

Rapid *E.coli*/Coliform Count Plates
Certificate #051801

AOAC® Official Method of AnalysisSM

Aerobic Count, Coliform Count Plates
AOAC 986.33, AOAC 989.10

High-Sensitivity Coliform Count Plates
AOAC 996.02

Aerobic Count Plates
AOAC 990.12

Coliform Count, *E.coli*/Coliform Count Plates
AOAC 991.14

Yeast and Mold Count Plates
AOAC 997.02

Rapid Coliform Count Plates
AOAC 2000.15

E.coli/Coliform Count Plates
AOAC 998.08

Enterobacteriaceae Count Plates
AOAC 2003.01

Salmonella Express System
AOAC 2014.01

Staph Express System
AOAC 2003.07, AOAC 2003.08, AOAC 2003.11

Rapid Yeast and Mold Count Plates
AOAC 2014.05

Rapid Aerobic Count Plates
AOAC 2015.13

AFNOR or MicroVal Certification (following ISO 16140-2)

Aerobic Count Plates
3M 01/01-09/89²

Rapid Aerobic Count Plates
3M 01/17-11/16²

Coliform Count Plates
01/02-09/89 A², 3M 01/02-09/89 B², 3M 01/02-09/89 C²

Select *E.coli* Count Plates
3M 01/08-06/01²

Rapid Coliform Count Plates (14 hour result)
3M 01/05-03/97 A²

Rapid Coliform Count Plates (24 hour result)
3M 01/05-03/97 B²

Enterobacteriaceae Count Plates
3M 01/06-09/97²

High-Sensitivity Coliform Count Plates
3M 01/07-03/99²

Staph Express System
3M 01/09-04/03 A², 3M 01/09-04/03 B²

Rapid Yeast and Mold Count Plates
3M 01/13-07/14²

Lactic Acid Bacteria Count Plates
3M 01/19-11/17²

Rapid *E. coli*/Coliform Count Plates
2017LR76

United States Industry Recognition

US FDA (United States Food and Drug Administration) Agricultural Marketing Service

Milk	Aerobic Count Plates Coliform Count Plates High-Sensitivity Coliform Count Plates Rapid Aerobic Count Plates 3M™ Petrifilm™ Plate Reader	<i>FDA Evaluation of Milk Laboratories</i> , 2017 Revision: https://www.fda.gov/media/115265/download
------	--	---

USDA FSIS (Food Safety and Inspection Service)

Examination of fresh, refrigerated and frozen prepared meat, poultry and pasteurized egg products	Aerobic Count Plates <i>E.coli</i> /Coliform Count Plates <i>Enterobacteriaceae</i> Count Plates	<i>Microbiology Laboratory Guidebook</i> , Chapter 3.01, Quantitative Analysis of Bacteria in Foods as Sanitary Indicators. January 20, 2011
---	--	--



3M™ Petrifilm™ Plate Reader Advanced


Coming soon!

To be the first to know, please visit us at [link
www.engage.3m.com/R3](http://www.engage.3m.com/R3)



Agenda

I. Evolution of Quantitative Microbial Methods



II. Evolution of Pathogen Detection Technologies



Foodborne pathogens

Pathogen: A bacterium, virus, or other microorganism that can cause disease.

- In most cases these microorganisms should be absent in a food or water sample
- Typically used as a qualitative assessment (presence/absence) to determine the microbiological safety
- For certain foodborne pathogens and particular food segments, a quantitative assessment may be necessary (Example: *Staphylococcus aureus*, *Bacillus cereus* and *Campylobacter*)

Common foodborne pathogens



Salmonella




Listeria monocytogenes



Campylobacter spp.



E. coli




Cronobacter

Other foodborne pathogens

- *Vibrio* spp.
- *Shigella*
- *Bacillus cereus*
- *Clostridium botulinum*
- Norovirus

Pathogen testing




To prevent foodborne illness

- Pathogen testing is an important element of food safety plans
- Increased trend in pathogen testing for incoming raw materials and food production environment


Measures to prevent contamination of end product

- Preventing the entrance of contaminated foods
- Monitoring the presence of indicators or pathogenic microorganisms in the food production environment
- Final product testing as part of the release criteria

Pathogen detection: A very brief timeline



Pathogen detection overview



Note: Rapid methods (molecular and immunoassay) give presumptive positive results

Methods of pathogen detection

Culture Method

Time consuming

Labor intensive

Look for specific metabolic traits

- Resistant to a particular antibiotic
- Ability to grow in the presence of a particular chemicals (salt, bile etc.)
- Ability to utilize a particular chemical or nutrient

Rely on SELECTIVE media to determine the presence of pathogens




Rapid

Results in less time

Reduce cost of operation with faster turnaround time

Immuno-based assays

- ELISA
- Immunochromatographic or dipstick assays



DNA-based assays

- PCR
- Isothermal DNA amplification



The 3M™ Molecular Detection System



Three main elements of 3M™ Molecular Detection System

The 3M™ Molecular Detection Assays use Loop-Mediated Isothermal Amplification (LAMP) to amplify specific regions of the DNA with high sensitivity and specificity.

1

An optimized enrichment protocol that will allow the growth and recovery of a target pathogen.



2

Specific 3M Molecular Detection Assays which provide all elements for sample preparation and DNA detection.









Salmonella
Listeria monocytogenes
Listeria spp.
E. coli O157:H7
Cronobacter
Campylobacter
STEC (2)

3

The 3M Molecular Detection Instrument for amplifying the DNA and the software to detect DNA amplification.



Simplified for Productivity

<p>3M Method <i>Salmonella</i></p> <p>37°C ±1°C or 41.5°C ±1°C 18-30 h</p> 	<p>3M Method <i>E. coli</i> O157 (including H7) STEC Gene Screen</p> <p>41.5°C ±1°C 10-24 h</p> 	<p>3M Method <i>Listeria</i></p> <p>37°C ±1°C 24-32 h</p> 	<p>3M Method <i>Listeria monocytogenes</i></p> <p>37°C ±1°C 24-32 h*</p>  <p>*Raw dairy matrices require 40 h</p>	<p>3M Method <i>Cronobacter</i></p> <p>37°C ±1°C 18-24 h</p> 	<p>3M Method <i>Campylobacter</i></p> <p>41.5°C ±1°C 22-28 h</p> 
--	---	--	--	--	--

One protocol for all
pathogen targets



Transfer 20 µL enriched sample to lysis tube.



Heat 15 minutes – 100°C ±1°C.
Cool 5 minutes on block at room temperature.



Transfer 20 µL lysate to reagent tubes
containing lyophilized pellet.



Place tubes in instrument.
Start run.
Amplification & detection in 15-75 minutes.
Automated & color-coded real-time results.

Dairy Processing




● *Listeria*




● *L. monocytogenes*



● *Salmonella*



● *Cronobacter*



USDA-FSIS MLG and FDA BAM Recognition



USDA FSIS Chooses 3M

USDA Food Safety and Inspection Service Chooses 3M for *Salmonella* and *Listeria monocytogenes* Testing.

3M Food Safety is honored to announce that following rigorous performance evaluations against other commercially available rapid methods, the United States Department of Agriculture's Food Safety and Inspection Service (USDA FSIS) has chosen [3M™ Molecular Detection System](#) as the primary method to be used for the detection of *Salmonella* and *Listeria monocytogenes*; two major pathogenic organisms continually threatening food production and processing.

[USDA FSIS MLG method 4.10: *Salmonella*](#)

[USDA FSIS MLG Method 8.11: *Listeria monocytogenes*](#)



Equivalent Testing Methodologies for *Listeria* species and *Listeria monocytogenes* in Environmental Samples

FDA has determined that the following methods are “scientifically valid” and “at least equivalent to the method of analysis in § 112.152(a) in accuracy, precision, and sensitivity”¹ in detecting *Listeria* species and *L. monocytogenes*. The method of analysis in § 112.152(a) is "[Testing Methodology for *Listeria* species or *L. monocytogenes* in Environmental Samples](#)" (October 2015, Version 1).

1. AOAC Official Method 2013.10. VIDAS® UP *Listeria* (LPT), *Listeria* in select foods and environmental surfaces.
2. AOAC Official Method 2016.07. 3M Molecular Detection Assay (MDA) 2, *Listeria* in select foods and environmental surfaces.
3. AOAC Official Method 996.14. Assurance polyclonal enzyme immunoassay method (EIA), *Listeria* in select foods and environmental surfaces.

Method Certifications & Validations

Recognized by organizations and government agencies around the world

International Recognition

AOAC® Performance Tested MethodSM

MDA2 – <i>Salmonella</i> ••Certificate #091501	MDA2 – <i>Listeria</i> ••Certificate #111501	MDA2 – <i>L. monocytogenes</i> ••Certificate #081501
MDA2 – <i>Cronobacter</i> ••Certificate #101703	MDA2 – <i>Campylobacter</i> ••Certificate #111803	MDA2 – STEC Gene Screen ••Certificate #071902 and 071903

AOAC® Official Method of AnalysisSM

MDA2 – <i>Salmonella</i> ••AOAC 2016.01	MDA2 – <i>Listeria</i> ••AOAC 2016.07	MDA2 – <i>L. monocytogenes</i> ••AOAC 2016.08	MDA2 – <i>E. coli</i> O157 ••AOAC 2017.01
MDA2 – <i>Cronobacter</i> ••AOAC 2018.01	<i>Salmonella</i> ••AOAC 2013.09	<i>L. monocytogenes</i> ••AOAC 2014.07	<i>Listeria</i> ••AOAC 2014.06

NF VALIDATION certificate granted by AFNOR Certification

MDA2 – <i>Salmonella</i> ••3M 01/16-11/16	MDA2 – <i>Listeria</i> ••3M 01/14-05/16	MDA2 – <i>L. monocytogenes</i> ••3M 01/15-09/16	
MDA2 – <i>Cronobacter</i> ••3M 01/20-03/18	MDA2 – <i>E. coli</i> O157 ••3M 01/18-05/17	<i>E. coli</i> O157 ••3M 01/12-03/13	<i>Salmonella</i> ••3M 01/11-11/12

USDA

MDA2 – <i>Salmonella</i> **MLG 4.10	MDA2 – <i>L. monocytogenes</i> **MLG 8.11
--	--

US FDA (Equivalent technology)

MDA2 – <i>Salmonella</i> ••AOAC 2016.01	MDA2 – <i>Listeria</i> ••AOAC 2016.07
--	--

Recognition by Country

Australia

Department of Agriculture/Australian Quarantine and Inspection Service
Approved Methods

MDA2 – <i>Salmonella</i> ••AOAC 2016.01	MDA2 – <i>Listeria</i> ••AOAC 2016.07
MDA2 – <i>L. monocytogenes</i> ••AOAC 2016.08	

Brazil

Ministry of Agriculture, Livestock and Supply (MAPA) Official Method

MDA2 – <i>Salmonella</i>	MDA2 – <i>L. monocytogenes</i>
--------------------------	--------------------------------

Central America

(Costa Rica, El Salvador, Guatemala, Honduras, Nicaragua and Panama)

Technical Regulation Annex to Resolution 402-2018

MDA2 – <i>Salmonella</i> ••AOAC 2016.01	MDA2 – <i>L. monocytogenes</i> ••AOAC 2016.08
MDA2 – <i>Cronobacter</i> ••AOAC 2018.01	MDA2 – <i>E. coli</i> O157 (including H7) ••AOAC 2017.01

Canada

Heath Canada Compendium of Analytical Methods

MDA2 – <i>Salmonella</i> ••MFLP-100	MDA2 – <i>Listeria</i> ••MFLP-101
<i>Salmonella</i> ••MFLP-06	<i>E. coli</i> O157 (including H7) ••MFLP-73
<i>Listeria</i> ••MFLP-05	<i>L. monocytogenes</i> ••MFLP-72



Thank you